

VacoSlides™

Integral Double Agar & Vacuum Test Tube

Automatic, one-step evacuation of low volumes (1/2 cc) for dipping & inoculating

Cross results to prevent errors

Start incubation at the point of collection


Easy to use

- quick, automatic operation
- minimal training required
- no opening containers or fumbling with accessories

Safe and convenient

- hermetically sealed — eliminates spills and clean-ups
- protects you and your staff from contamination during handling and transfer
- urine odors minimized from sample collection to the lab

Urine collection combined with culture media




Available with CHROMagar/MacConkey/CLED or others

User features


- designed for use with BD Vacutainer, Greiner, or Agit Urine Cups
- no need to open VacoSlides test tube or Urine Cup
- automatic, just insert the test tube

Lab testing features

- fast visual inspection of negative results (75% of the tests)
- fast visual bacteria identification and colony density
- easy opening for continued sensitivity test, if necessary



To Order:
Agit Technologies (Israel) www.Agit.co.il



**Association
for Public
Health Services**
For a Better Quality of Life