

From urine collection straight to the lab

VacoLoop™

Integral Inoculating Loop & Vacuum Test Tube

Extract micro samples for agar streaking
Quickly, cleanly, effectively

Easy to use

- easy sample collection and retrieval in any setting

Saves money

- a fraction of the cost of other sample collection systems

Reduces false negatives & false positives

- reliable and accurate results every time

Spill-Proof

- integral system design minimizes spills and clean-ups

Biohazard safe


- protects you and your staff from contamination and infection

Urine Odors Minimized

- minimal odors from sample collection to the lab

VacoLoop Features

- Designed for use with BD-Vacutainer, Greiner, and Agit Urine Cups
- Accurately deposits liquid, covering only the VacoLoop eye
- Immediately and safely transfers self-contained sample to the lab
- Automatically evacuates liquid directly from the urine cup
- Streaks volumes of 1 µL or 10 µL automatically and accurately every time
- CE and AMAR (Israel Medical Product Certification) Certified; ISO 13485 Certified facility




Patient ID barcode

STUDIO SHUKI KOOK
March 11 - 6:20 BROCOIE

Ordering Information	VL10	Vacoloop 10 µL, sterile
	VL01	Vacoloop 1 µL, sterile

To Order:
Agit Technologies (Israel) www.Agit.co.il



Association for Public Health Services
For a Better Quality of Life